

# EYE DRAPE

SURGIWEAR has developed this drape after feedback from surgeons. It offers many advantages:

- Level 3 protection as per European Union Norms.
- Operating area has 2 layer protection.
- No fluid can pass through the drape in operating area.
- Minimum transfer of fabric fibers into the wound.
- No chemicals present in the fabric of the drape.
- Incision drape is made of breathable lo-Drape.

Use of this drape will result to:

- No or less post operative swelling of wound.
- Negligible post operative fever
- Supple wound area.
- Minimum post operative pain.
- Lesser incidence of post operative infections.
- Lesser morbidity and faster recovery of patient.
- **Ultimate cost saving by way of lesser use of medicines.**

**Uses:** For all kind of surgery on Eye

## e EYE DRAPE D700

**Material :** 22-30 gsm  
SMS polypropylene fabric water repellent  
Sheet Size : 100 x 70 cm  
Incise Area : 7 x 9 cm  
Eye lid holders : ( 6 no)  
Drainage Pouch of  
Transparent PE film: (1 no)  
Poly sheet lining 60 x 60 cm  
Sterilized by EO



## e EYE DRAPE ECO E711

**Material :** 25 gsm  
SMS Polypropylene fabric water repellent  
Sheet Size : 100 x 70 cm  
Incise Area : 7 x 9 cm  
Eye Lid Holders : ( 6 no )  
Drainage Pouch of Transparent PE film : ( 1 no )  
Sterilized by EO



## e EYE DRAPE LARGE D702


**Material :** - 65 - 77 gsm  
Spun Lace fabric water repellent  
Sheet Size: - 120 x 160 cm  
Incise Area : 17 x 7 cm  
Drainage Pouch of Transparent PE film ( 1 no )  
Eye Lid Holders ( 6 no )  
Sterilized by EO

## e EYE DRAPE BLUE D703

**Material :** 40-50 gsm  
SMS polypropylene fabric water repellent  
Sheet Size : 100 x 75 cm  
Incise Area : 7 x 9 cm  
Hole in Drape : 3 x 5 cm  
Eye lid holders : ( 6 no)  
Drainage Pouch of Transparent PE film : ( 1 no)  
Poly sheet lining 60 x 60 cm  
Sterilized by EO

## e EYE DRAPE NW BLUE LARGE D704

**Material :** 40 -50 gsm  
SMS Polypropylene fabric water repellent  
Sheet Size : 100 x 160 cm  
Incise Area : 7 x 9 cm  
Hole in Drape : 3 x 5 cm  
Eye Lid Holders : ( 6 no )  
Drainage Pouch of  
Transparent PE film: (1 no )  
Sterilized by EO



## e EYE DRAPE POLY BLUE LARGE D705


**Material:** Blue PE film 25 micron  
Sheet Size: - 100 x 120 cm  
Incise Area : - 7 x 9 cm  
Hole in Drape : - 3 x 5 cm  
Drainage Pouch of Transparent PE film ( 1 no )

be all free from disease & sufferings

**SURGIWEAR®**  
EN ISO 13485: 2012 COMPANY



Eye Drape NW Blue Large D704



Eye Drape NW Blue

Eye Lid Holders ( 6 no )  
Sterilized by EO

### *m* MINOR DRAPE PLAIN D600

**Material :** PE Hygiene film 25 micron  
Sheet Size : 60 x 60 cm  
Incise Area : 6 x 8 cm  
Eye lid holders : ( 6 no )  
Sterilized by EO

### *e* EYE DRAPE POLY D601

**Material :** PE Hygiene film 25 micron  
Sheet Size : 60 x 60 cm  
Incise Area : 6 x 8 cm  
Eye Lid Holders : ( 6 no )  
Drainage Funnel of  
Transparent PE film: (1 no )  
Sterilized by EO



Eye Drape Poly

### *e* EYE PATCH D602

**Material :** 25 micron PE film Adhesive coated with  
40 - 60 gsm adhesive.  
Sheet Size: 13 x 15 cm  
Drainage Pouch of Transparent PE film ( 1 no )  
Eye Lid Holders ( 6 no )  
Sterilized by EO

**SURGIWEAR®**  
EN ISO 13485: 2012 COMPANY

[www.surgiwear.co.in](http://www.surgiwear.co.in)

Made in India by:

**G. Surgiwear Limited**

Rasoolpur Jahanganj, Shahjahanpur-242001

Tel: +91 5842-223818 Fax: +91 5842-222190

e.mail: [surgiwear@hotmail.com](mailto:surgiwear@hotmail.com)

Delhi: 65454565, Mumbai: 66994258, Ahmedabad: 27545218