


# SHUNT PARTS & ACCESSORIES

## S "SLIT N SPRING" VALVE



High Pressure	SH211
Medium Pressure	SH212
Low Pressure	SH213

### "Slit N Spring" Valve with Large Flushing reservoir



High Pressure	SH221
Medium Pressure	SH222
Low Pressure	SH223

### Anti-Siphon / Anti- Gravity Valve Chhabra "Z" Flow Valve



Two Ball System	SH111
Three Ball System	SH112
Four Ball System	SH113

## u UNI-SHUNT

90 cm long radio-opaque tube OD 2.5mm & ID 1.3 mm. Proximal end bears 32-40 holes. The tip is closed by a tantalum impregnated plug. Distal end is closed with silicone plug and bears four vertical slits. The pack also contains two straight connectors. Supplied ready to use in double peel open packs, sterilized by EO.

High Pressure	SH401
Medium Pressure	SH402
Low Pressure	SH403

## d DISTAL PART SH001



90 cm long radio-opaque tube OD 2.5mm & ID 1.3 mm. Distal end is closed with silicone plug and bears four vertical slits. Supplied ready to use in double peel open packs, sterilized by EO.

## v VENTRICULAR PART

### Ventricular Part SH002



15 cm long radio-opaque tube OD 2.5mm & ID 1.3 mm. The end bears 32-40 holes. The tip is closed by a tantalum impregnated plug. Supplied ready to use in double peel open packs, sterilized by EO.

be all free from disease & sufferings

**SURGIWEAR®**  
EN ISO 13485: 2012 COMPANY



#### **Ventricular Part Right Angle 5 cm SH002R5**

15 cm long radio-opaque tube OD 2.5mm & ID 1.3 mm. 5 cm long ventricular end. The end bears 32-40 holes. The tip is closed by a tantalum impregnated plug. Supplied ready to use in double peel open packs, sterilized by EO.

#### **Ventricular Part Right Angle 7 cm SH002R7**

15 cm long radio-opaque tube OD 2.5mm & ID 1.3 mm. 7 cm long ventricular end. The end bears 32-40 holes. The tip is closed by a tantalum impregnated plug. Supplied ready to use in double peel open packs, sterilized by EO.

#### **Ventricular Part SH002B**

30 cm long radio-opaque tube OD 2.5mm & ID 1.3 mm (75F). The end bears 32-40 holes. The tip is closed by a tantalum impregnated plug. Supplied ready to use in double peel open packs, sterilized by EO.

#### **Ventricular Part SH002C**

30 cm long radio-opaque tube 10.5F. The end bears 32-40 holes. The tip is closed by a tantalum impregnated plug. Supplied ready to use in double peel open packs, sterilized by EO.

#### **LUMBAR PART SH009**

30 cm long radio-opaque OD 1.5 mm ID 0.7 mm . Proximal end bears multiple holes. The tube is marked at 5 cm, 10 cm and 15 cm distance from proximal perforated tip. It can pass through 14 gauge Touhey needle. The catheter is supported by a teflon coated guide wire with flexible tip. Supplied ready to use in double peel open packs sterilized by EO.

#### **a ANTE CHAMBER SH008**



For connection with Chhabra shunt. For CSF sampling or injection of medications into ventricles. Can be punctured by 27 gauge needle at an angle 10-15 times.

#### **S SUTURE COLLAR SH007**



Used to anchor the tubing with thread. It is made of silicone.

#### **C CONNECTOR**

**Y connector with extra long ventricular part SH006**



Made of Kynar. Supplied with Ventricular part 30 cm long radio-opaque tube OD 2.5mm & ID 1.3 mm (75F). The end bears 32-40 holes. The tip is closed by a tantalum impregnated plug. Supplied ready to use in double peel open packs, sterilized by EO.

#### **Straight Connector SH003**



Either made of PEEK or Kynar. Supplied ready to use in peel open packs, sterilized by EO.

#### **Right Angle Connector SH004**



Made of Kynar. Supplied ready to use in peel open packs, sterilized by EO.

## Lumbar Catheter Connector SH005



Made of Peek & stainless steel. Supplied ready to use in peel open packs, sterilized by EO.

## C CEREBRAL RESERVOIR



Small 12 mm dia

SH011

Large 18 mm dia

SH012

Reservoir Chamber is made of silicone and base is made of poly propylene. Supplied with one Ventricular part 15 cm long radio-opaque OD 2.5 mm ID 1.3 mm closed at the end within black silicone tantalum impregnated plug. The end has 32-40 holes supplied with stylet. Supplied ready to use in double peel open packs sterilized by EO.

## i INSTRUMENTS



Sub Cutaneous Catheter Passer Small (35 cm) INS01, Large (60 cm) INS02



Ventricular Catheter Passer INS03



Peritoneal Split Trocar INS04

### PRODUCT INFORMATION DISCLAIMER

G. Surgiwear Limited has exercised reasonable care in the choice of materials and manufacture of this product. G. Surgiwear Limited excludes all warranties, whether expressed or implied by operation of law or otherwise, including, but not limited to any implied warranties of merchantability or fitness for a particular purpose. G. Surgiwear Limited shall not be liable for any incidental or consequential loss, damage or expense, directly or indirectly arising from use of this product. G. Surgiwear Limited neither assumes nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this product.

be all free from disease & sufferings

**SURGIWEAR®**  
EN ISO 13485: 2012 COMPANY

**SURGIWEAR®**  
EN ISO 13485: 2012 COMPANY

[www.surgiwear.co.in](http://www.surgiwear.co.in)

Made in India by:

**G. Surgiwear Limited**

Rasoolpur Jahanganj, Shahjahanpur-242001

Tel: +91 5842-223818 Fax: +91 5842-222190

e.mail: [surgiwear@hotmail.com](mailto:surgiwear@hotmail.com)

Delhi: 65454565, Mumbai: 66994258, Ahmedabad: 27545218